

ACTIVALL™

Dynamic standing frame



ATTENTION!


Manufacturer takes responsibility for proper usage of the device only when the product was purchased from qualified Akces-MED Ltd. employee or at specialized medical store!

Akces-MED Ltd. reserve the right to introduce technical and trading changes in the content of the instruction without warning.



Product marked with the CE sign is Class I medical device, of rule 1, consistent with the requirements of Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017 on medical devices.

Akces-MED Ltd. company introduced quality management system ISO 13485 regarding production, sales and service.



Contents

- Introduction.....4
- 1. Purpose of user manual.....5
- 2. Indications and contraindications for use.....5
- 3. General characteristic.....7
- 4. Size table.....8
- 5. How to assembly the product.....9
- 6. How to use the product.....11
 - 6.1 Chest support.....11
 - 6.2 Pelvic support.....12
 - 6.3 Pelvic harness.....13
 - 6.4 Castors with friction brake.....14
- 7. Learning to walk.....15
- 8. The use of the product.....16
 - 8.1 Preparation of the user.....16
 - 8.2 Preparation of the product.....17
 - 8.3 Upright position.....18
- 9. Safety rules.....19
- 10. Terms of warranty.....21
- 11. Maintenance and cleaning.....21
- 12. Labelling.....25

Introduction

As far this dynamic standing frame is concerned, it is highly beneficial to the process of daily rehabilitation which takes place not only in hospitals and rehabilitation wards but also in user's house. This device enables the user to stand and walk. Moreover, children are able to take first steps and adults can either keep or regain fitness.

Both light and stable frame is helpful in independent relocation in this device with slight involvement of the user. The correct adjustment of braces makes the user feel secure as well as comfortable during performing activities. Furthermore, this product may be used by users in different age groups due to a vast array of sizes.

ACTIVALL™ plays a significant role because of the fact that user has the possibility of gaining new experience and, more to the point, being independent as far as possible.

Remember to read this user manual carefully. You will find not only useful information but also guidelines regarding safe and effective use of the product.

If you have any questions or remarks, do not hesitate to contact with us:

Medical consultant: mobile +48 516 183 883

Akces-MED Ltd.: ph. +48 17 864 04 77

e-mail: export@akces-med.com

ATTENTION!

Before using the product, read this user manual which is essential for proper and safe usage of the product.

Remember that complying with guidelines is highly important. It assures safety and extends durability as well as aesthetics of the dproduct.

1. Purpose of user manual

ATTENTION!

The present user manual contains basic information which is essential for: preparation before use, proper usage, maintenance, cleaning and guarantee conditions.

User manual should be kept in easily accessible place. It is designed for people who look after the disabled, and for both doctors as well as physiotherapists who operate the product.

Before using the ACTIVALL™ dynamic standing frame user is obliged to familiarize himself / herself with user manual.

2. Indications and contraindications for use

The standing frames are intended for people with orthopedic, neurological and cardiological diseases and for people after surgery who are unable to adopt and maintain a standing position on their own as a result of prolonged immobilization. The device is an irreplaceable aid in the rehabilitation process, both in terms of movement and social integration. It is great help in everyday activities such as rehabilitation but also during playing or studying. The selection of the appropriate size of the device depends on medical indications and anatomical parameters, regardless of the age of the user.

It is recommended for users in the following cases:

- cerebral palsy – various types and forms,
- brain and spinal cord injuries with lower limbs paralysis and weakness, balance and motor coordination disorders with different etiology,
- demyelinating disease (e.g. multiple sclerosis – SM),
- meningocele,
- muscles diseases with weakness and paralysis (dystrophy, myasthenia gravis, myositis),
- genetic disorders,
- degeneration and neurological diseases,
- other condition connected with muscles weakness, paralysis and musculoskeletal disorders,
- after head and spinal cord injuries and spine trauma with paralysis and paresis.

2. Indications and contraindications for use

Contraindications:

While standing upright, the user who is in a standing position should be closely monitored. The therapy should be discontinued in case of such side effects as: turning pale, sweating or fainting, joint strain and the resulting deformities of the body. Further proceeding should be consulted with the medical specialist or a physiotherapist who prescribed the treatment.

The critical condition of the user is an absolute contraindication to standing upright.

Relative contraindications for upright standing are:

- arterial hypertension at rest above 200/110 mmHg, pulmonary hypertension, cardiac arrhythmias, unstable asthma, unstable metabolic disorders, moderate heart valve disease,
- any condition, such as a healing fracture or severe osteoporosis, that makes it impossible to bear the user's own bodyweight,
- bone fracture with incomplete union, unless approved by a doctor,
- postoperative restrictions on joint loading,
- pain while standing in the standing frame (caused by subluxation or dislocation of the hip joint, overstretching of soft tissues or excessive pressure). Subluxation or dislocation of the hip joint in the absence of pain is not a contraindication.

People with significant flexion in the hip or knee joints, as well as with contractures in the ankle joint, may not be able to use the standing frame. Make sure that the position of the user in the standing frame does not cause excessive pressure or stretching of affected areas of the body. In such cases, standing upright may be contraindicated until a greater range of motion is achieved in the contracted joints.

In case of people with impaired cardiovascular system (arterial hypertension at rest above 200/110 mmHg, pulmonary hypertension, cardiac arrhythmias, moderate heart valve disease); respiratory system (unstable asthma) and unstable metabolic disorders, it is necessary to monitor the user while standing in order to prevent serious complications.

3. General characteristic

ACTIVALL™

- ADJUSTABLE CHEST SUPPORT
 - ADJUSTABLE PELVIC SUPPORT
 - PELVIC HARNESS
 - HAND GRIP
-
- FRAME WITH CASTORS EQUIPPED WITH BRAKES
 - SEPARATOR OF THE LOWER LIMBS* OR STABILIZER OF THE LOWER LIMBS**
 - FRONT CASTORS WITH BRAKES - 75 mm
 - REAR CASTORS WITH FRICTIONAL BRAKES - 75 mm




- * STABILIZER OF THE LOWER LIMBS - it converts ACTIVALL™ into a fully functional static version of the stander.
- ** SEPARATOR OF THE LOWER LIMBS - it prevents the legs from crossing during the use of the device's dynamic function

INFORMATION!

Available accessories can be found on the website akces-med.com in the product section.

4. Size table



			Size				
Symbol	Measurement	Unit	1	2	3	4	5
W	indicative user height *	[cm]	70 - 100	100 - 130	110 - 150	110 - 170	140 - 200
A / A1	pelvic / chest support width	[cm]	26	28	31	37	42
C	pelvic support height	[cm]	44 - 60	66 - 82	77 - 102	83 - 112	98 - 120
C1	chest support height	[cm]	52 - 70	77 - 95	89 - 116	92 - 126	108 - 143
B / B1	pelvic / chest support depth *	[cm]	22	25	27	29	35
	pelvic / chest support perimeter	[cm]	87	97	106	124	144
	max. user weight	[kg]	25	40	60	80	95

Weight and dimensions		Unit	1	2	3	4	5
width		[cm]	59	60	63	66	74
length		[cm]	73	84	84	91	98
height		[cm]	54	79	90	96	112
weight		[kg]	13,5	16,5	18	20,5	23

Attention: Measurement tolerance range is +/- 1 cm / ½ in.
Attention: Approximate values, depending on the product configuration.
* Measured without pelvic harness.
** The value cannot be main factor for selection of product size.

5. How to assembly the product

INFORMATION!

Before starting the assembly of the product, make sure that the surface where the product is being mounted is flat and does not pose a risk of the stander tipping over during assembly.

It is recommended to activate the brakes after removing the base.


ATTACHING THE COLUMNS TO THE BASE OF THE STANDER

The columns should be placed on the base at the hole locations and screwed in using the screws provided with the columns. This step does not apply to size 1 of the stander, where the columns are permanently attached to the base.


INSTALLING THE PELVIC SUPPORT

Insert the supports of the hip clamp into the first frame of the base. After placing the clamp and adjusting it to the desired height, tighten the knobs securely.

ATTACHING THE COLUMNS TO THE BASE OF THE STANDER



INSTALLING THE PELVIC SUPPORT




5. How to assembly the product

INFORMATION!

INSTALLING THE CHEST SUPPORT

Wsunąć wsporniki obejmą piersiowej w drugą ramę podstawy. Po założeniu obejm i ustawieniu jej na wybranej wysokości, należy mocno dokręcić pokrętła.

INSTALLING THE CHEST SUPPORT



ATTENTION!

Before using the product, it is essential to check the correctness of its assembly.

Ensure that the screws securing the stander columns are properly tightened and that the knobs are securely fastened.

Additionally, check if the hip and chest clamps do not shift under pressure.

6. How to use the product

ATTENTION!

The initial adjustment of the stander to the user should be carried out solely by properly trained personnel, and further adjustments should be made by trained individuals. Proper adjustment of the product is essential for the safety and ensures the correct course of the user's rehabilitation.

Before using the stander, make sure to lock the wheel brakes. Unlocked wheels may cause the product to shift during user placement, potentially leading to a fall and resulting in injury.

Before placing the user in the stander, the supports should be opened to ensure easy access to the product. It is important not to rest the body weight on the open clamp protections. These are not intended to assist during entry or exit from the stander. Overloading them may cause deformation of the parts, malfunction of the snap lock, and loss of product stability, potentially leading to injury.

Once the user is fully upright, close the clamps, ensuring they are securely locked.

6.1 Chest support

INFORMATION!

The adjustment of the chest support is made by loosening the knobs on the stander columns, setting the desired height of the support, tightening the knobs firmly, and ensuring they are correctly tightened.

Attention should be paid to avoid exceeding the double lines on the clamp supports.



ATTENTION!

The height of the chest support depends on the user's control over the chest. If the control is good, the user can maintain an upright torso position without significant effort, feeling secure and comfortable - then the support can be positioned slightly lower.

However, the walking process should begin with the support placed under the armpits, followed by gradual lowering as needed.

6. How to use the product

6.2 Pelvic support

PELVIC SUPPORT HEIGHT ADJUSTMENT



INFORMATION!

The adjustment of the hip support is done by loosening the knobs on the stander columns, setting the desired height of the support, tightening the knobs firmly, and ensuring they are correctly tightened.

Attention should be paid to avoid exceeding the double lines on the clamp supports.

PELVIC SUPPORT DEPTH ADJUSTMENT



INFORMATION!

To adjust the depth of the support, loosen and remove the knobs located beneath the support, set the desired position (so that the holes on the frames align), and then screw the knobs back in, ensuring they are correctly and securely tightened.



6. How to use the product

6.2 Pelvic support

ATTENTION!

The pelvic support should be positioned at the level of the greater trochanters of the lower limbs.

Po ustawieniu wysokości obejmę biodrowej należy dopasować stopień jej wysunięcia przód - tył względem bioder. Od dokładności dopasowania obejmę biodrowej zależy prawidłowa pozycja tułowia i miednicy.

6.3 Pelvic harness

INFORMATION!

The pelvic harness secures the user against falling.

It is fastened using buckles on the pelvic support.

Before use, ensure that the buckles do not unfasten under the force applied to the harness.



ATTENTION!

The pelvic harness should not be treated as a seat. It is intended solely as a safety measure to prevent the user from falling.

INFORMATION!

The upright position is both safe and comfortable. In this position, the user adapts to the vertical posture, the sense of balance adjusts, and the internal organs and circulatory system adapt to new conditions and increased load.

In this position, the user can remain as long as needed while performing exercises, small tasks, etc.

The product remains stable even in the event of the user fainting.

6. How to use the product

6.4 Castors with friction brake

INFORMATION!

The use of rear wheels with a resistance brake allows for the regulation of the user's movement speed in the stander.

This is particularly important when starting to learn to walk in the stander.

The knob is used to adjust the pressure applied to the brake pad, which slows down the wheel's rotation.

When taking the first steps in the ACTIVALL™, it is recommended to use the resistance brake to make the device move more slowly. Once the user gains confidence, the brake resistance can be reduced.



TIGHTENING THE KNOB

- INCREASES THE PRESSURE OF THE BRAKE PAD ON THE WHEEL
- SLOWS DOWN THE SPEED OF THE STANDER



LOOSENING THE KNOB

- REDUCES THE PRESSURE OF THE BRAKE PAD ON THE WHEEL
- INCREASES THE SPEED OF THE STANDER



7. Learning to walk

The process of learning to walk in a dynamic stander is a gradual rehabilitation process. The dynamic stander is a device that supports correct posture and enables a step-by-step approach to walking.

Preparing the user for the walking process

A doctor or physiotherapist evaluates the user's physical and neurological condition to tailor the device to their needs, considering posture, joint range of motion, and muscle strength.

The device is then adjusted according to the user's individual parameters, such as height, weight, and the need for support in different body parts.

Positioning the user

The initial phase involves learning to maintain a standing posture in the stander. The user, supported by caregivers, becomes accustomed to the new position while simultaneously working on balance and full body stabilization. The time spent standing in the stander is gradually extended based on the user's capabilities.

Activating leg muscles

Once the user is able to maintain an upright position for an extended period, exercises to activate the leg muscles are introduced. Gradually, more advanced exercises are introduced to prepare the user for attempting to take steps.

Preparing for walking

The user, gradually supported by the device, starts taking their first steps. A crucial part of the therapy is practicing balance and coordination during walking attempts, learning the rhythm of walking, and synchronizing the movements of the lower limbs with the rest of the body.

The dynamic stander helps maintain balance and provides cushioning in case of weakness or fainting while walking. At this stage, as the user takes steps and gains control over their body, both the distance covered and the duration of the exercises are gradually increased.

It is important that the walking learning process is tailored to the user's individual needs and abilities.

After completing exercises with the stander, the user, with the help of a caregiver, should safely return to a wheelchair, chair, or bed.

ATTENTION!

The user should never leave the stander independently. The caregiver is responsible for supporting the user from the side or behind.

Particular attention should be given to locking the wheels of the wheelchair or stabilizing the chair or other equipment before moving it backward during the user's transfer. The back movement of the seat could lead to the user falling, resulting in injury.

It is prohibited to move the user in the stander across sloped, uneven surfaces or areas with obstacles (e.g., curbs, steps), as this could damage the device.

8. The use of the product

ATTENTION!

Before using the ACTIVALL™ dynamic standing frame, read this user manual.

Before each use, prepare the user and the device!

8.1 Preparation of the user

Before using the product, it is essential to make sure that:

- user's outfit is light and comfortable sporty clothes are advisable,
- shoes should be comfortable and, moreover, they should stabilize ankle,
- due to hygienic reason and to avoid abrasion it is advisable for the user to always wear clothes when using the device, skin must not be naked especially in places of regular contact with upholstery,
- drains, bags, catheters etc. if applied, have to be secured to function properly and not to disturb during using the product,
- if it is necessary for a user to wear orthopedic supply (collar, corset, orthosis, shoes) if it is recommended by doctor.

ATTENTION!

Outfit should protect user's skin from direct contact with upholstery as well as other parts of equipment.

Only doctor can make a decision concerning using any kind of orthopedic supply!

8. The use of the product

8.2 Preparation of the product

Before use it is necessary to:

- assemble the standing frame,
- place the standing frame on an even and flat surface, away from sources of fire and high temperature,
- make sure that all the elements are undamaged, properly mounted and stable,
- push the brakes,
- adjust the standing frame to the user's size and needs,
- mount accessories (if applicable),
- place the user in the standing frame,
- fasten the elements which stabilize the user,
- start uprighting.

ATTENTION!

Before using the stander, ensure that it has been correctly assembled, and that the knobs, screws, and other fastenings are properly tightened. Check for any loosening of the fastenings under applied load.

Before using the stander, verify that all straps and safety devices have been properly adjusted and fastened.

The stander will only function correctly if it is tailored to the user's dimensions and needs. The primary rule when adjusting the stander is to provide the user with maximum safety. Do not force any body position changes, as this may cause discomfort. Improper use of the stander can lead to the reinforcement of incorrect posture patterns.

The initial adjustment of the stander to the user should only be performed by appropriately trained personnel, and continued by trained individuals. Proper fitting of the device is essential for safety and ensures the correct course of the user's rehabilitation.

8. The use of the product

8.3 Upright position

The length of time for upright position process is different for each user. It is recommended to extend this time during the initial period (adaption period for the change of position).

Either doctor or physiotherapist can make decision concerning the length of upright position process. Upright position process has to be always performed under the supervision of trained person!

It is recommended to divert user's attention during upright position process by: manipulation games, educational games, touch stimulation and others.

ATTENTION!

It is highly important to adjust the device to the user before using it!

Before starting the standing process, ensure that all supports and the pelvic harness have been properly adjusted and fastened, and that the fastenings of adjustable elements are securely tightened.

It is forbidden to leave the user unattended!

Only trained staff can adjust standing frame to the user!

As far as first upright standing is concerned, it takes place in rehabilitation, neurosurgery or orthopedics wards. Moreover, it has to be performed together with doctor and physiotherapist due to the fact that medical problems may appear during first treatments such as fainting, drop in blood pressure, dizziness etc. Upright standing is continued in home together with trained staff (rehabilitation technician, district nurse, trained family members).

9. Safety rules

ATTENTION!

The dynamic standing frame ACTIVALL™ has to be used in compliance with guidelines in this use manual.

Before each use, prepare the user and the device!

It is highly important to consult either doctor or physiotherapist in order to make right decision concerning the choice of ACTIVALL™ and its size. It is also essential to familiarize with the guidelines described in user manual before using the device.

1. Before use it is necessary to:
 - make sure that there are no damages or different faults which could endanger user's life,
 - check if all parts are working and are correctly mounted and, moreover, check if they are not broken,
 - check if all screws, knobs and all glued, sewed or tapped parts are tightened,
 - check the condition of the upholstery for any signs of damage,
 - check if brakes are operating correctly,
 - check if buckles are operating correctly,
 - adjust the device according to user's measurements and needs e.g. the adjustment of supports, harness and accessories (if any),
 - make sure that all movable parts and knobs for adjustments are stable and locked,
 - secure the user with available stabilizing aids of the standing frame. It is highly important to secure patient because it safeguards him/her from falling out.
2. It is forbidden to use standing frame for different purposes that it is intended.
3. It is forbidden to use standing frame when it is damaged or any part is missing.
4. It is forbidden to leave the user unattended. Device should be used only when there is someone who is trained and knowledgeable.
5. It is forbidden to either go up or down the stairs. It doesn't matter if the user is placed in the device or not. It is forbidden to leave the user and the device on unstable ground, even if the brake is engaged.
6. It is important to engage brake during getting into and getting out the standing frame.
7. Make sure that all limbs are not close to spokes or other adjustable elements.
8. The standing frame should be stored and used away from sources of fire or high temperature. Moreover, standing frame shouldn't be stored in rooms with high humidity.
9. Secure standing frame during transportation in order to avoid any possible damages. Manufacturer doesn't take responsibility for damages which appeared due to improper preparation of the device during transportation.
10. Before using the standing frame, remove all paper and plastic packaging. The packaging is not a toy! Keep away from children and animals.

9. Safety rules

ATTENTION!

Incorrect product adjustment can lead to such side effects as fainting, body deformities, joint overload, etc. If the abovementioned symptoms occur, stop using the product immediately, take the user out of the product and lay them down. If the undesirable symptoms persist, contact a physician immediately. If the symptoms have disappeared, inform the attending physician about this fact in order to make a decision about the further rehabilitation process.

Serious incidents involving this device should be reported both to manufacturer and to the competent authority in your country.

The device contains small parts which children may swallow.

Do not use a product that shows any signs of damage, destruction or wear. Only a fully functional product, used in accordance with the instructions for use, is safe for the user. If the product is damaged, discontinue use and contact us with the manufacturer.

It is forbidden to make any changes concerning device's structure due to the fact that it is dangerous and it results in losing warranty.

The weight of user should not be bigger than maximal weight described on page 8 where you will find information regarding this issue.

10. Terms of warranty

1. Akces-MED Ltd. sells products free from manufacturing defects or faults.
2. Manufacturer gives the purchaser 24-month guarantee from the date of purchase. All defects revealed during warranty period will be repaired free of charge.
3. All defects which appear after warranty period can be repaired only when the user will pay for it.
4. The stander is covered by a 36-month warranty on the territory of Canada – the details of the rights resulting from the extended warranty are regulated by the warranty card.
5. Manufacturer claims that the product can be used for 48-months. User can use this product after this period only when he/she will send the device to manufacturer so as to service it. Next, manufacturer determines again how long the product can be used and, moreover, the date of the next service will be determined.
6. If product is damaged, please contact our dealer or Akces-MED Ltd. customer service: (ph. +48 17 864 04 77, e-mail: export@akces-med.com).

11. Maintenance and cleaning

The dynamic standing frame ACTIVALL™ which is a medical device, should always be kept clean and should be used in accordance with the manufacturer's instructions. The standing frame should be kept in a clean, dry and ventilated room. Maintenance and cleaning must be done at least once a month. Periodic cleaning and maintenance of the device is essential to keep its proper functioning and safety.

Recommendations for cleaning the upholstery and other elements of the device:

- wash with warm water using a soft cloth and mild cleaners,
- make sure the cloth is not too wet, as this may enlarge the stain and damage the upholstery or foam inside the upholstery,
- remember to remove excess moisture from the upholstery before the next cleaning,
- the upholstery must not get wet during cleaning,
- you must not soak the upholstery,
- after cleaning the stains, the upholstery should be rinsed with lukewarm water and dried with a clean cloth,
- do not use hard brushes for cleaning the upholstery,
- do not clean the frame and the upholstery with strong chemicals, e.g. solvents or chlorine-based caustics,
- do not clean, replace any elements or maintain the standing frame during use,
- metal and plastic elements should be cleaned by hand with suitable cleaners, after having read the recommendations of the manufacturer of the given cleaner,

11. Maintenance and cleaning

- it is possible to disinfect metal, plastic, wooden and upholstery elements of the product with a generally available disinfectant. In order to disinfect the selected element, spray its surface with a disinfectant generally available on the market and wipe with a dry, clean cloth. After disinfection, wait until the element is completely dry, especially in the case of upholstery. Before using the disinfectant, first check in an invisible place whether the product will not discolor or damage the disinfected surface.

Before using it is necessary to:

- clean wheels,
- check if wheels are correctly mounted, check if brake stops turning wheel,
- check if standing frame construction is free of breaks or any kinds of deformities,
- check if all movable parts are working and all screws, nuts, knobs etc. are correctly tightened,
- check if belts are not frayed; check if buckles and clamps work correctly.

ATTENTION!

Remember that all elements and upholstery has to be always dry before using.

The manufacturer of the dynamic standing frame, Akces-MED Ltd., is not responsible for any damages caused by using unsuitable cleaners.

Manufacturer takes responsibility only for hidden faults or damages which appeared despite the fact that user complied with the guidelines.

In case of any damages, it is essential to contact the person who sold you the product.



11. Maintenance and cleaning

TECHNICAL INSPECTION

Technical inspections should be performed regularly. Before servicing, the product should be thoroughly cleaned to reveal any hidden damage.

FREQUENCY	WHAT SHOULD BE CHECKED?	WHAT SHOULD BE DONE	REMARKS
BEFORE EACH USE			
	General condition of the product.		
	Check:		
	<ul style="list-style-type: none">the entire product for any visible damage,cracks or mechanical defects,if the product complete,if all bolts, screws, knobs, nuts and other elements of the construction are complete and tightened.		
	Wheels and brakes.		
	Check:		
	<ul style="list-style-type: none">if the wheels are properly mounted and if they are not damaged,if the brake is blocking the spinning wheel,if the wheels are not dirty. If necessary, clean the wheels.		
	Fastening of the construction elements.		
	Check all the mechanisms of the product and its moving parts. Check in particular:		
	<ul style="list-style-type: none">if the product is properly mounted,if the wheels are properly mounted,if the elements securing the user are properly mounted and if they work well,if all the moving parts are complete,if all the moving parts work properly, do not jam,if the adjustment mechanisms of the following elements work properly (if applicable).		
	Upholstery elements.		
	Check:		
	<ul style="list-style-type: none">if there is any damage, tears, rips, loose seams or deformations on the upholstery elements,if the pelvic harness and straps are not frayed or twisted and if the buckles work properly,if the palvic harness do not unfasten under pressure or when pulled.		
	In case of any damage or defects, discontinue the use of the product and contact your Supplier / Distributor or the Manufacturer.		
ONCE A WEEK			
	Product cleaning.		
	If the frame and upholstery get dirty (dusty, muddy, etc.) during everyday use, wipe them with a damp cloth without detergents. Remove dirt (e.g. hair, food remains, etc.) from the moving parts of the product.		
	The product should be cleaned when necessary, but at least once a week. Do not use any chlorine based cleaners or methylated spirit.		
	Before using the product, check if the upholstery is dry.		

























11. Maintenance and cleaning

TECHNICAL INSPECTION

Technical inspections should be performed regularly. Before servicing, the product should be thoroughly cleaned to reveal any hidden damage.


FREQUENCY	WHAT SHOULD BE CHECKED?	WHAT SHOULD BE DONE	REMARKS
ONCE A MONTH			
	Wheel axles, quick connectors (if applicable).		
	Check: <ul style="list-style-type: none">if the wheels are rolling and if they are not loose. To maintain the wheels, apply a small amount of a preservative-lubricant between the axle and the wheel hub. Check if you can easily remove and attach the wheels. Clean the wheels when dirty.		
	Product construction connectors (welding, threading, soldering).		
	Check: <ul style="list-style-type: none">the condition of the product frame in the places of joining (welds, threads, soldering).		
	Tires/innertubes (if applicable).		
	Check: <ul style="list-style-type: none">tire wear condition. Check for visible damage . Check tire pressure (if applicable). If the tire pressure is too low, inflate them to the correct level.		
	In case of any damage or defects, discontinue the use of the product and contact your Supplier / Distributor or the Manufacturer.		
ONCE EVERY 6 MONTHS			
	Construction moving parts.		
	Check: <ul style="list-style-type: none">if the moving parts of the construction have not been damaged. Maintain them using aproper cleaning agent.		
	Product labelling.		
	Check: <ul style="list-style-type: none">if the marking and labels are readable.		
	In case of any damage or defects, discontinue the use of the product and contact your Supplier / Distributor or the Manufacturer.		
AFTER THE LIFE TIME OF THE PRODUCT FINISHES			
	Detailed product check.		
	The manufacturer specifies the life time of the product, which is given in the user manual. After the given period, the user may use the product only if it is sent back to the manufacturer to be checked first. After that, the manufacturer will re-define the life time of the product and the date of the next check.		

12. Labelling

	Declaration that the product meets the essential requirements		Dry vertically
	Serial number		Do not chloriante
	Please read the manual instructions before use		Do not iron
	Attention!		Do not dry clean
	Manufacturer		Do not bleach
	Production date		Carefully, delicate content
	The product does not contain latex		It is forbidden to either go up or down the stairs. It doesn't matter if the user is placed in the device or not
	Store in a dry place		Prohibition to leave the user without care during use of the product
	Protect from sunlight		Prohibition of use without fastened safety belts
	Storage temperature limitation		You need to read the additional information
	Medical device		You need to read the additional information
	Unique Device Identification		You need to use the tools indicated in the assembly description

The product data plate contains basic product information such as name, maximum load, manufacturer's data, date of manufacture, and marking in accordance with the requirements for medical devices.

MANUFACTURER



Jasionka 955 B
36-002 Jasionka
ph. +48 17 864 04 7
e-mail: export@akces-med.com
BDO: 000031619

Authorised Dealer