

PARAPION™

Static standing frame

We know very well how important it is for the user to adopt the standing position, that is why we have created a product called PARAPION™ which is ideally suited for this purpose. PARAPION™ allows you for full and safe upright standing. Stable construction, ergonomic knee supports, wide pelvic and chest belts are the basic elements ensuring safety and comfort during the therapy. The user, who is placed in this device, is able to perform rehabilitation treatment as well as simple domestic or professional activities. The standing frame is equipped with four noiseless castors which enable easy relocation of the user.



Easy operation

24⁺

24-month warranty



Stable construction

Colors

Frame color



Size 1



Size 2



Size 3




70 - 130 cm

Indicative user height



35 - 45 kg

Max. user weight



For **more details** about the product, please visit our website.

This is a medical device.
Use it in accordance with the user manual or the label.

INDICATIONS

Cerebral palsy (CP) – various types and forms | Brain and spinal cord injuries of various etiologies with limb paresis, severe balance disturbances, and coordination impairments | Demyelinating diseases (multiple sclerosis – MS, amyotrophic lateral sclerosis – ALS, and others) | Meningocele | Muscle diseases with paresis and paralysis (including dystrophy, myasthenia gravis, polymyositis) | Genetic syndromes with limb paresis | Degenerative and metabolic diseases of the nervous system | Degenerative and rheumatic joint diseases | Other diseases with limb paresis, paralysis, and musculoskeletal impairments | After injuries, following traumatic brain injuries or spinal cord injuries with limb paresis and paralysis, particularly with quadriplegia (as part of the rehabilitation process or as permanent orthopedic support)

CONTRAINDICATIONS

Observation: Monitor user during standing. Stop therapy if paleness, sweating, fainting, joint overload, or body deformations occur; consult with a doctor or physiotherapist | **Absolute Contraindication:** Critical condition of the user. | **Relative Contraindications:** Resting hypertension > 200/110 mmHg, pulmonary hypertension, heart rhythm disorders, unstable asthma, unstable metabolic issues, moderate heart valve disease. | Conditions like healing fractures or severe osteoporosis that prevent weight-bearing | Incomplete fracture healing (unless approved by a doctor) | Postoperative joint load restrictions | Pain from dislocation, excessive tissue stretch, or pressure; dislocation without pain is not contraindicated | Significant hip/knee flexion or ankle contractures; ensure no excessive pressure or stretch on affected areas.